



# MYTHIC™ 5VET PRO

Fully automated 5DIFF hematology analyzer for veterinary use.



## Technical specifications

|                           |  |
|---------------------------|--|
| Throughput                | 60 samples per hour  |
| Sample volume             | 20 µL (whole blood and pre-diluted mode)   |
| Pack of reagents          | Reagents Diluent, Optical Fluid, Lysing Reagent, Cleaner reagents recognized with RFID code  |
| Technology                | <ul style="list-style-type: none"> <li>• Laser Light Multi-Dimensional Cell</li> <li>• Classification (WBC DIFF)</li> <li>• Impedance (WBC, RBC, PLT)</li> <li>• Spectrophotometry (HGB)</li> </ul>        |
| 34 parameters             | WBC, LYM% and #, MON% and #, NEU% and #, EOS% and #, BASO% and #, RBC, HGB, HCT, MCV, MCH, MCHC, RDW_CV, RDW_SD, PLT, MPV, PDW, PCT, P_LCC, P_LCR, RET%, RET_ABS, IRF, NRBC% and #, LIC% and #, ALY% and # |
| Histograms & Scattergrams | RBC histogram, PLT histogram, Scattergrams for WBC differentiation, Scattergrams for RET   |
| Dimensions                | 270 × 508 × 412 mm (W x D x H)   |
| Weight                    | 26 kg  |
| Power                     | <ul style="list-style-type: none"> <li>• 100 V to 240 V</li> <li>• 50 to 60 Hz</li> <li>• Less than 250 VA</li> </ul>  |
| Printers                  | <ul style="list-style-type: none"> <li>• Built-in thermal printer</li> <li>• External printer USB port</li> </ul>  |
| User's interface          | <ul style="list-style-type: none"> <li>• 10.4" touchscreen</li> <li>• info about keyboard to be deleted</li> </ul>   |
| Output                    | <ul style="list-style-type: none"> <li>• RS232</li> <li>• Ethernet (UDP/IP, TCP/IP protocol)</li> <li>• 2 USB port</li> <li>• Centronics port</li> <li>• Barcode scanner port</li> </ul>                   |

## Software specifications

|                       |  |
|-----------------------|--|
| Graphic interface     | Intuitive menu with big icons  |
| Available species     | Dog, Cat, Horse, Pig, Cow, Buffalo, Rabbit, Monkey, Rat, Mouse, Sheep, Goat, Camel, 3 User-defined   |
| Multilanguage         | English, Chinese, Spanish, French, Russian, Portuguese, Czech, Turkish, Polish, Italian  |
| Calibration           | Automatic or manual  |
| Quality control       | <ul style="list-style-type: none"> <li>• 6 batches: 100 run/batch</li> <li>• Table and Levey Jennings graphs</li> </ul>  |
| Archives              | <ul style="list-style-type: none"> <li>• Internal memory: 200 000 patients</li> </ul>  |
| Traceability          | <ul style="list-style-type: none"> <li>• Patient identification: name, patient ID, sample ID</li> <li>• Daily log</li> <li>• Reagents monitoring: lot number, date, volume</li> <li>• Service log</li> </ul> |
| Flagging              | <ul style="list-style-type: none"> <li>• Pathologic flags</li> <li>• Lab limits</li> <li>• Reagent alerts</li> <li>• Instrument alerts</li> </ul>  |
| Reagents<br>HMYET-010 | Mythic 5 VET PRO Diluent<br>Mythic 5 VET PRO Optical Fluid<br>Mythic 5 VET PRO Lysing Reagent<br>Mythic 5 VET PRO Cleaner<br>MYTHIC RETIC PRO  |



Professional Veterinary use Hematology system for the complex 5DIFF blood count with Reticulocytes determination

- i** Orphée reserves own rights to technical and software modification.  
In case of printer choice please contact your local Orphée's representative.  
Number of tests per pack depends on the analyzed species and daily workflow.  
Mass storage depends on capacity of USB mass storage device.



Orphée S.A. 19, chemin du Champ-des-Filles  
CH-1228 Plan-Les-Ouates, Geneva, Switzerland

phone: +41 22 884 90 90  
fax: +41 22 884 90 99

www.orphee-medical.com  
contact@orphee-medical.com